



25mm Filter
Nylon 0.45µ

25mm Filter
Nylon 0.22µ

25mm Filter
PTFE 0.45µ

25mm Filter
PTFE 0.22µ

25mm Filter
RC 0.45µ

Color Coded Filters

Easy to Identify*, Low Volume Syringe Filters for HPLC, CE or Dissolution Testing

*Membrane Type and Pore Size Imprinted on Filter

Low Cost

Available in:

- Aqueous Compatible: Nylon
- Chemically Inert: PTFE
- Low Protein Binding: RC

Supplied with a Certificate of Analysis
Each Filter Tested for Membrane Integrity

Specifications

Membranes: Nylon, PTFE and Regenerated Cellulose (RC)
Filter Size: 25mm
Retention Volume: 90µl
Maximum Operating Pressure: 110psi
Maximum Operating Temperature: 100°C
 Device is made from Virgin Medical Grade Polypropylene

For more information visit www.mtc-usa.com

Cat. No.	Description	Pore Size	Qty.	Free Samples*
58045-N25-A	Nylon, Hydrophilic	.45µ	50	
58045-N25-C	Nylon, Hydrophilic	.45µ	100	
58045-N25-K	Nylon, Hydrophilic	.45µ	1,000	
58022-N25-A	Nylon, Hydrophilic	.22µ	50	
58022-N25-C	Nylon, Hydrophilic	.22µ	100	
58022-N25-K	Nylon, Hydrophilic	.22µ	1,000	
58045-P25-A	PTFE, Chemically Inert	.45µ	50	
58045-P25-C	PTFE, Chemically Inert	.45µ	100	
58045-P25-K	PTFE, Chemically Inert	.45µ	1,000	
58022-P25-A	PTFE, Chemically Inert	.22µ	50	
58022-P25-C	PTFE, Chemically Inert	.22µ	100	
58022-P25-K	PTFE, Chemically Inert	.22µ	1,000	
58045-R25-A	RC, Low Protein Binding	.45µ	50	
58045-R25-C	RC, Low Protein Binding	.45µ	100	
58045-R25-K	RC, Low Protein Binding	.45µ	1,000	

Q. How do I get samples of these filters to try?

A. Three easy ways...

1. **Internet.** Visit our website and fill out our "Sample Request Form" under Customer Support menu.
2. **Phone.** Give us a call, our customer service department will gladly walk you through the fast process.
3. **Fax.** Send us a fax on your company letter head, including your shipping address, fax and phone and we will fax back a confirmation.

* Free to qualified/approved scientists.

25mm 1-05



MTC • MicroSolv Technology Corporation

1 Industrial Way West, Bldg E, Eatontown, NJ 07724 USA

Ph. 1-732-578-1777

Fx. 1-732-578-9777

www.mtc-usa.com

Email: info@mtc-usa.com